

## Interi System Product List

The Interi System Kit includes everything needed to initiate a patient treatment and send the patient home using Interi. The appropriate Kit should be selected and made available in the OR at the time of surgery.

Interi Kit contains all of the system components required for surgical placement:

- Manifold
- Manifold Connector
- Therapy Units
- Interi System Instructions for Use Healthcare Professional

Interi Kit also contains items patients will need to use Interi at home:

- Extra Therapy Units
- Therapy Unit carrying case
- Go-home bag to transport the items upon discharge
- Interi System Instructions for Patient Use

Interi Therapy Unit 4-Packs contain individually packaged Therapy Units which can be on the shelf at the surgeon's office, clinic or hospital. This allows flexibility to dispense Therapy Units to patients on an as-needed basis should they require additional units to manage their fluid removal needs.

SKU	Name	150mL Therapy Units per Kit	300mL Therapy Units per Kit	500mL Therapy Units per Kit	Total Volume Capacity
K300450	Interi Kit 3-450	3			450mL
K300600	Interi Kit 3-600	2	1		600mL
K300900	Interi Kit 3-900		3		900mL
K301100	Interi Kit 3-1100		2	1	1100mL
K301500	Interi Kit 3-1500			3	1500mL
P400600	Interi Therapy Unit 4-Pack - 150mL	4			600mL
P401200	Interi Therapy Unit 4-Pack - 300mL		4		1200mL
P402000	Interi Therapy Unit 4-pack - 500mL			4	2000mL